

	Configurations — must be ordered with the station											Options — may be added in the field		Add at any time
	 Cryogen-free	 Magnetic field	 Pneumatic vibration isolation	 Stand	 12.5:1 microscope	 High vacuum	 Load-lock [†]	 Very low temperature operation	 Backside optical access	 360° sample stage rotation	 ±5° sample stage rotation	 High temperature operation	 Low temperature operation	 Ring magnet sample holder kit
TTPX*			✓	✓	✓				✓			to 675 K	to 3.2 K	Up to ~0.19 T
CPX			✓	included	✓	✓	✓	to 1.6 K			■	to 675 K	to 1.9 K	Up to ~0.19 T
CPX-VF		±2.5 T vertical	✓	included	✓	✓	✓				■		to 1.9 K	
CPX-HF		±1 T horizontal	✓	included		✓					■		to 1.9 K	
EMPX-H2		±0.6 T horizontal		included						✓			to 3.2 K	
FWPX			included	included							■		to 3.5 K	
CRX-6.5K	■			included	✓							to 675 K		Up to ~0.19 T
CRX-4K	■			included	✓							to 675 K		Up to ~0.19 T
CRX-VF	■	±2.5 T vertical		included	✓	✓	✓							
CRX-EM-HF	■	±0.6 T horizontal		included	✓					✓				

■ = feature inherent to system ✓ = feature available blank cell = feature not available

*The TTPX is available in a fully specified, preconfigured tabletop model with four-arm triaxial configuration (PS-100).

[†]As a field-installable option to the load-lock, a suitcase for controlled sample transport is also available.